

NexTemp[®] Ultra

Disposable Thermometer

- **1 FULL MINUTE SIGNAL RETENTION TIME**

Longer time to read temperature.

- **THROW AWAY THE GERMS**

Reduces the risk of cross-infection or reinfection.

- **CONVENIENT**

Easily transportable, unbreakable, inexpensive.

- **CLINICALLY ACCURATE**

$\pm 0.2^{\circ}\text{F}$ / $\pm 0.1^{\circ}\text{C}$ (within the critical range of 98.0°F to 102.0°F / 37.0°C to 39.0°C).

- **LONG SHELF LIFE**

5 years from the date of manufacture.

- **SAFE**

Non-toxic, environmentally friendly, latex-free.


- **ECONOMICAL**

No probe covers, batteries, calibration checks, capital equipment.

- **NO SPECIAL STORAGE CONDITIONS**

Resets automatically if exposed to elevated temperatures during shipment or storage. Best if stored in cool dry place (86°F / 30°C or below).

DIRECTIONS FOR USE – See box or package insert for details.



Last black dot shows 98.6°F

Last black dot shows 37.2°C

HOW TO READ THE THERMOMETER



Ordering Information

ULTRA	$^{\circ}\text{F}$	8000-15	1/2 Case	10 Boxes of 100 (1,000 pcs.)
	$^{\circ}\text{F}$	8000-20	1 Case	20 Boxes of 100 (2,000 pcs.)
	$^{\circ}\text{C}$	8200-15	1/2 Case	10 Boxes of 100 (1,000 pcs.)
	$^{\circ}\text{C}$	8200-20	1 Case	20 Boxes of 100 (2,000 pcs.)

Individually Wrapped
Disposable ULTRA
Clinical Thermometers

KwikTemp Medical

UEI: HAC5GD8G44A6 • CAGE Code: 8MC97
246 North New Hope Road, Gastonia, NC 28054

Phone: 980-253-0949 • Contact KwikTemp • www.kwiktempmedical.com

